


MediSystem Pro1 Multifrequency Body Composition Analyser Nutritional and Health Evaluation

Incluyes 5 Software Screens:

1. Personal Data



MediSystem Pro1

File Data System Help

Personal Data | Body Composition | Parameters | Fitness | Reports

State: Active

Personal Data | Factors

Name: Andrew Surname: JENKINS File num: 1

Address: 46 Fairchild Crescent Town: Midtown

City: Southampton Post Code: SO9 - 5YT


Phone 1: 16733992 Phone 2: Mobile: Fax:

Date of Birth: 11/12/1959 Identif: E-mail: landy@hotmail.com

Additional Info 1: Additional Info 2:

Comments:

2. Body Composition – Clic on “Measure Impedance”



MediSystem Pro1

File Data System Help

Personal Data | Body Composition | Parameters | Fitness | Reports

Anthropometrical Data

Sex: Male

Population: Adult

Build: Medium

Age: 39 years

Weight: 92 Kgr

Height: 179 cm

4 Electrodes Bioimpedance

Read Impedance

1khz	10khz	25khz	50khz	75khz	100khz	150khz
658	634	621	612	590	570	550

Measurement Procedure

1. Click 'Measure Impedance'
2. Click 'Measure'
3. Wait for end of measurement
4. Click 'Save Measurement'
5. Click 'Read Impedance'
6. Select the File added at point 4.

Measure Impedance

3. Parameters: Circumferences, Cardiac Function & Biochemical

MediSystem Pro1
File Data System Help

Personal Data | **Body Composition** | **Parameters** | Fitness | Reports

Perimeters

Chest cm
 Thigh cm
 Waist cm
 Hip cm

Cardiac Function

Systolic pressure mmhg
 Diastolic pressure mmhg
 Pulse beats/min

Biochemical

Total Cholesterol mg%
 HDL Cholesterol mg%
 LDL Cholesterol mg%
 Triglicerides mg%
 Iron mcg%
 Glucose mg%
 Uric acid mg%
 Urea
 Creatinin mg%
 GOT - AST ui

GPT - ALT ui
 GAMMA - GT ui
 Hematin millon/ul
 Leukocyte millon/ul
 Hemoglobin g/dl
 Hematocrit %

MSPROWEBSUPPORT-ENG - Microsoft Word

4. Fitness: Fitness Planner & Other Tests

MediSystem Pro1
File Data System Help

Personal Data | Body Composition | Parameters | **Fitness** | Reports

Stamina

Pres ups
 Abdominals

Flexibility

Sit and Reach Cms
 Sapo Cms

Muscle strength

Jump Cms

Cardiac resistance


Rest beats/min
 Immediatly after test beats/min
 1 minute after test beats/min

Occupation


Fitness Planner

(Minutes)	Monday	Tuesda	Wednes	Thursda	Friday	Satur
Badminton	60			60		
Basketball						
Canoeing						
Cycling 8,8k						
Cycling 15,1						
Dancing						
Hockey (gra						
Golf		180			180	
Gymnastics						
Horseriding						
Jogging 9mi						
Judo, Karate						
Running 6 m						
Running						
Ski-ing cross						
Ski-ing down						
Squash						
Swimming sl						


5. Reports – Calculates Results Immediately




Report Configuration: Choose the Content and Cover Page



Example Report – Easy to Read with Colour Graphs



Preview Your Report Before Printing



Export Results and Parameters of All Files to Open in Excel Spreadsheet



	D	E	F	G	H	I	J	
1	Weight	Height	Complexion	Population	Weight	Ideal	Fat	Le
2								
3		92	179	Medium	Adult	92 69,14- 74,56	39,5	52
4		69	178	Medium	Adult	69 68,58- 74,02		12

MediSystem Pro1 is very easy to Install and Use

1. Install the Software



2. Contact the Analyser to the RS232 (serial) port or USB port of computer

